

PALACOS® fast R+G

The new fast setting cement for experienced surgeons

PALACOS® fast R+G

- Fast setting cement solution based on the proven bone cement PALACOS® R+G
- Fast setting cement with total working time less than 5 minutes (at room temperature)
- For permanent fixation of joint endoprotheses, especially acetabular cup and small joints
- New member of the PALACOS® product family – PALACOS® fast R+G



Product	Composition	Content	REF
PALACOS® fast R+G	Fast setting, high viscosity bone cement with Gentamicin	1x40	66056768

Simply order PALACOS® fast R+G from Heraeus.

Australia

orders:

Toll Free: 1800 226 521

phone: +61 (02) 8422 6100

email: sales.hkaus@heraeus.com

www.heraeus-medical.com

06720

PALACOS® fast R+G

Save time: no waiting time – short setting time

The new member in the PALACOS® product family – PALACOS® fast R+G

PALACOS® fast R+G is a fast setting, high viscosity bone cement, containing Gentamicin. It is based on the well known and proven bone cement PALACOS® R+G and offers reduced working times.


PALACOS® fast R+G is sterilized via Ethylenoxide gassing which leads to a high reproducibility of handling and mechanical properties.¹

Optimal application time for experienced surgeons

PALACOS® fast R+G offers an application timeframe that meets the need for the application in small joints and acetabular cups in hip arthroplasty.

PALACOS® fast R+G was developed mainly for experienced surgeons and can be used to fix total or partial joint endoprostheses.

Working time comparison of PALACOS® fast R+G and PALACOS® R+G



Test conditions: Not pre-chilled bone cement, manually mixed at room temperature (23°C)

¹ Kühn, KD. Bone Cements. Up-to-Date Comparison of Physical and Chemical Properties of Commercial Materials. Berlin: Springer, 2000; 33.